

Self-contained System for Mitigation of Contaminated Aerosol Sources of SARS-CoV-2

by

Adithya Shyamala Pandian

A Thesis Presented in Partial Fulfillment
of the Requirements for the Degree
Master of Science

Approved April 2021 by the
Graduate Supervisory Committee:

Erica Forzani, Chair
Bhavesh Patel
Xiaojun Xian

ARIZONA STATE UNIVERSITY

May 2021

ABSTRACT

Contaminated aerosols and micro droplets are easily generated by infected hosts through sneezing, coughing, speaking and breathing¹⁻³ and harm humans' health and the global economy. While most of the efforts are usually targeted towards protecting individuals from getting infected,⁴ eliminating transmissions from infection sources is also important to prevent disease transmission. Supportive therapies for Severe Acute Respiratory Syndrome Coronavirus 2 (SARS CoV-2) pneumonia such as oxygen supplementation, nebulizers and non-invasive mechanical ventilation all carry an increased risk for viral transmission via aerosol to healthcare workers.⁵⁻⁹ In this work, I study the efficacy of five methods for self-containing aerosols emitted from infected subjects undergoing nebulization therapies with a diverse spectrum on Non-Invasive Positive Pressure Ventilator (NIPPV) with oxygen delivery therapies. The work includes five study cases: *Case I*: Use of a Full-Face Mask with biofilter in bilevel positive airway pressure device (BiPAP) therapy, *Case II*: Use of surgical mask in High Flow Nasal Cannula (HFNC) therapy, *Case III*: Use of a modified silicone disposable mask in a HFNC therapy, *Case IV*: Use of a modified silicone disposable mask with a regular nebulizer and normal breathing, *Case V*: Use of a mitigation box with biofilter in a BiPAP. We demonstrate that while cases I, III and IV showed efficacies of 98-100%; cases II and V, which are the most commonly used, resulted with significantly lower efficacies of 10-24% to mitigate the dispersion of nebulization aerosols. Therefore, implementing cases I, III and IV in health care facilities may help battle the contaminations and infections via aerosol transmission during a pandemic.

ACKNOWLEDGMENTS

I would like to thank my advisor, Dr. Erica Forzani for her supervision and support throughout the work. I would have not been able to complete my work without her help. I would like to thank Dr. Bhavesh Patel for helping in facilitating the experiments and providing valuable guidance throughout the work. I would also like to thank Dr. Xiaojun Xian for taking out the time to serve on my committee.

I am thankful to my colleagues at Biosensors Lab, Bio-design Center at ASU for their support and guidance. I am deeply thankful to Dr. Kelly McKay for extending his support by participating as the subject in the experiment. I would also like to thank Dr. Amelia Lowell for helping in conducting the experiment. I am happy to give my thanks to Mayo Clinic, Phoenix to provide the necessary infrastructure and equipment to carry out the experiments

Finally, I am very grateful to my dad for supporting me always. I would also want to thank my mom for her endless love and believing in me. I am always indebted to her a lot for all her sacrifices for me. I am also deeply thankful to my uncle for guiding me to be a better man.

TABLE OF CONTENTS

	Page
LIST OF FIGURES.....	v
CHAPTER	
1 INTRODUCTION	1
Research Motivation	1
Need for Better Mitigation System.....	1
2 NEBULIZATION OF A SARS-CoV-19 PATIENT	3
Expiremental Setup	3
Observations.....	4
3 MITIGATION STRATERGY	6
4 METHODS	8
Aerosol Sensing Instruments.....	8
Verification Test of the Aerosol Sensing Instruments	8
Experimental Rooms	10
Real and Simulated Human Subjects	12
Oxygen Delivery Methods	13
Aerosol Mitigation Methods	13
5 MITIGATION SYSTEM EVALUATION	14
Validation Test of the Aerosol Sensing Instruments	14
Assessment of Air Exchange Rate of the Rooms	15
Description of Study Cases	16
6 RESULTS	33
7 CONCLUSION	37
REFERENCES	38

LIST OF FIGURES

Figure		Page
1.	Nebulization Experimental Design	4
2.	Nebulization Aerosol Dispersion	5
3.	Different Cases of Combination of Oxygen Therapy and Mitigation System	7
4.	0.19 μ m Physiological Solution Particle Count 5 Minutes after Injection	9
5.	Experimental Testing Plan to Assess the Efficacy	12
6.	0.19 μ m Physiological Solution Particle Count.....	16
7.	Particle Count of BiPAP Therapy, with Full-Face Mask	19
8.	Particle Count of HFNC Therapy, with and without Surgical Mask	23
9.	Particle Count of HFNC Therapy, with and without Modified Silicon Mask.....	28
10.	Particle Count of Natural Breathing, with and without Modified Silicon Mask.....	33
11.	Particle Count of BiPAP Therapy, with Mitigation Box	38
12.	Efficacy of Aerosol Mitigation	44
13.	Modified Silicon Mask	47

CHAPTER 1

INTRODUCTION

1.1 Research Motivation

The coronavirus (COVID-19) pandemic has already infected over 45 million across the world and is responsible for more than 1.1 million deaths as of November 1st 2020.¹⁰ The disease is caused by the Severe Acute Respiratory Syndrome Coronavirus 2 (SARS-CoV-2).¹¹ The mortality rate is estimated to be nearly 5% for those between the ages of 45-64 years old and a striking 19.1% for those in the 75+ age group.¹² Clearly, SARS-CoV-2 is a massive risk given the widespread transmission of the virus.

1.2 Need for Better Mitigation System

To date, there is a plethora of evidence demonstrating that SARS-CoV-2 transmission occurs fundamentally from the spread of viral pathogens in the infected host's respiratory system to other susceptible hosts in contact with droplets, aerosols and fomites.^{3, 7, 13, 14} An aerosol is a suspension of fine particles (which can include viral pathogens) in an airborne liquid mist and can be transported through ventilation systems (e.g. AC) since it is not strongly affected by gravity.³ Since SARS-CoV-2's effective size is ~100 nm,^{11, 15} it can be encapsulated within aerosols from the respiratory system of an infected person, presenting a major health risk to the environment of any building where a SARS-CoV-2 infected person might reside or be present.^{1, 3} Ninety-nine percent of aerosols produced by humans, regardless of age, sex, weight and height are less than 10 μm .^{2, 3} This is concerning since the smaller the aerosol, the longer it takes to settle increasing the risk of inhaling the contaminated aerosols by other individuals.¹⁶ For example, an 0.5 μm aerosol takes 41 hours to settle.¹⁷ In case of SARS-CoV-2, the viruses can be viable on a surface for up to 3 days.¹⁸ Researchers have found SARS-CoV-2¹⁴ with virulent activity³ in collected aerosol particles from 0.2 μm to 10 μm , which is a serious concern for the spread of the disease through

air conditioning (AC) systems.^{7, 19} More recently CDC has recognized the importance of the transmission potential of SARS-CoV-2 via aerosols.²⁰ The groundbreaking evidences strongly indicates the need of aerosol mitigation to safeguard the public and public spaces of the populations.

CHAPTER 2

NEBULIZATION OF A SARS-CoV-19 PATIENT

2.1 Experimental Setup

In the present work, we target the mitigation of aerosol dispersion during nebulization treatments to provide a means of safe treatment for a respiratory disease such as COVID-19. This treatment causes a high risk of spreading pathogens due to the nebulization therapy generating aerosols with sizes less than $10 \mu\text{m}^2$, with emphasis on aerosol portions that do not impact the alveolar area, but remain in the dead space of the respiratory system (including nose and mouth) in contact with infectious areas and are exhaled into the environment⁹. The measurement was carried out in a Room of 17' x 13.6' x 9', equipped with a small bathroom of 4' x 6.25' x 9' in the back and an air ventilation rate, which is purposely set at ~ 20 air changes per hour⁻¹ to minimize transmission of the disease via aerosols^{17, 19}. The measurement was performed with a total of 3 aerosol sensors located at a height of 3 to 5 feet at 3 (DL-1), 6 (DL-2) and 9 feet (DL-3) away from the subject as shown in **Figure 1**. The sensors were able to detect aerosol particles equal to or larger than $0.19 \mu\text{m}$ (see **Figure 4**). The COVID-19 patient was assisted through an oxygen delivery therapy via a high flow nasal canula (HFNC) at 60 L/min. A nebulization therapy and an exercise therapy were delivered to the patient under the supervision of a respiratory therapist during the measurement. The nebulization included a total 3 ml solution containing 2.7 ml of Albuterol and 0.3 mL of saline physiological solution delivered through a Piezoelectric Based Nebulizer²². The subject did not wear any mask to mitigate the aerosol dispersion.




Figure 1. Nebulization Experimental Design. Experimental design for testing the nebulization’s aerosol dispersion in the patient’s room (air exchange rate = 20 h^{-1}).

2.2 Observations

Figure 2 illustrates the problem by showing aerosol particle counts / feet^3 transient over time inside a room with a COVID-19 patient during a nebulization therapy. As it can be observed, significant particle counts / feet^3 were detected above the baseline level of the room, typically at $\sim 7,400$ counts / feet^3 (for $0.2 \mu\text{m} - 5\mu\text{m}$ size count). During nebulization, the peak level was 365 times greater than the baseline level, and the particles clear out from the room 20 min after the completion of the nebulization. In addition, the proactive respiratory therapy executed by a therapist and patient’s breathing exercises under various conditions did not produce detectable aerosol-particle counts. Further, the action of flushing the bathroom toilet did not produce a particle count on the closest particle counter located 6 feet from the bathroom since the bathroom was small and

had its own ventilation system with air exchange rate of $\sim 20 \text{ h}^{-1}$. This seminal data indicates the nebulization therapies in COVID-19 patients should be targeted as a main source of potential contaminated aerosol in the environment.




Figure 2. Nebulization Aerosol Dispersion. Particle count concentration during HFNC treatment with 21% of Oxygen, 30 L/min of air flow rate and nebulization with 3 ml of medicine in saline physiological solution to a COVID-19 patient, respiratory exercise therapy and toilet flushing activity. The aerosol particle count profiles are simultaneously taken at 3 different positions of: 3 feet (—), 6 feet (---) and 9 feet (···) from the subject by an aerosol particle counter system (DL). The subject has no mitigation mask.

CHAPTER 3

MITIGATION STRATEGY

Aerosol spreading patterns are dependent on air exchange rates and ventilation streams of the room, the mechanism of nebulizers, and environmental (e.g. temperature and humidity) conditions. In fact, **Figure 2** indicates that the spatial distribution of the aerosol in the room is counterintuitive. Although the aerosol concentration at 3 feet was the highest, the peak concentration at 9 feet was 1/3 larger than the corresponding peak at 6 feet location. In order to solve the problem of highly variable and unpredictable aerosol patterns, we mitigate the dispersion by *DIRECTLY attacking the problem at the point of contaminated dispersion source and studying the effect of different mitigation systems that could potentially self-contain nebulizer aerosol dispersions.* The mitigation systems include 1- a biofilter which is a bacterial/viral filter with a humidity exchanger, 2- a silicone disposable mask (Breezing™) modified with a biofilter and a built-in fan, 3- a surgical mask, 4- a mitigation box with a biofilter. The application of the mitigation methods to different oxygen therapies rendered five study cases, as follows:




Figure 3. Different Cases of Combination of Oxygen Therapy and Mitigation Systems.

Case I: Use of a Full-Face Mask with biofilter in bilevel positive airway pressure device (BPAP) therapy (**Figure. 3-I**), Case II: Use of surgical mask in High Flow Nasal Cannula (HFNC) therapy (**Figure. 3-II**), Case III: Use of a modified silicone disposable mask in a HFNC therapy (**Figure. 3-III**), Case IV: Use of a modified silicone disposable mask with a regular nebulizer and normal breathing (**Figure. 3-IV**), Case V: Use of a mitigation box with biofilter with Non-Invasive Positive Pressure Ventilator (NIPPV) (**Figure. 3-V**).

CHAPTER 4

METHODS

4.1 Aerosol sensing instruments

The particles and aerosols generated during the nebulization procedures are assessed using the following commercial optical particle counters: MET ONE HHPC2+ from Beckman Coulter, and Dylos DC1100 Pro and Dylos DC1700 from Dylos Corporation, CA. MET ONE HHPC2+ particle counter provides reading for two ranges of particle sizes: 0.5 μm and 5.0 μm , which typical correspond to measuring particles $\geq 0.5 \mu\text{m}$ and $\geq 5.0 \mu\text{m}$, respectively. Dylos DC1100 Pro and Dylos DC1700 particle counters provide readings for particle sizes of 0.5 μm and 2.5 μm , which typical correspond to measuring particles $\geq 0.5 \mu\text{m}$ and $\geq 2.5 \mu\text{m}$, respectively. The sensors from Dylos DC1100 Pro and Dylos DC170 have identical sensing chambers but differ in the way they are powered. While Dylos DC1100 Pro is powered with a power adapter, Dylos DC1700 can be battery operated.


To decide on the best particle counters for this study, two additional particle counters were evaluated. These include: Dylos DC1100 from Dylos Corporation, CA with a capacity to detect to two ranges of particles: 1.0 μm and 5.0 μm (representing detection to $\geq 1.0 \mu\text{m}$ and $\geq 5.0 \mu\text{m}$, respectively), and Dusttrak DRX aerosol monitor from TSI Incorporated with a capacity to detect Aerosol concentration range 0.001 to 150 mg/m^3 .

4.2 Verification test of the aerosol sensing instruments

Evaluation of all above-mentioned particle counters were performed in a sensing chamber using pure polystyrene particles (Bang Laboratories, Inc.) of size 0.19 μm suspended in ethanol. The particles were aerosolized with the mechanical nebulizer and introduce to a sensing chamber with the aid of a small fan integrated in the chamber inlet port. Different particle aerosol concentrations were created by aerosolizing different volumes from 0.1 mL to 3.0 mL. The

sensitivity of the different particle counters to detect 0.19- μm polystyrene aerosol particles was evaluated and are shown in **Figure 4**. From all the particle counters evaluated, MET ONE HHPC2+ and Dylos DC1100 Pro demonstrated the capacity to detect 0.19 μm polystyrene particles with the highest sensitivity of all four counter, and therefore, they were the aerosol sensors of choice for the study.

For sake of simplification, we refer to MET ONE HHPC2+ particle counter as “MO sensor”, and to Dylos DC1100 and DC1700 as “DL sensor”.






Figure 4. 0.19 μm Physiological Solution Particle Count. (A) The Experimental setup with four types of particle counter – i. MetOne, ii. Dylos 0.5 μm , iii. Dylos 0.1 μm and iv. Dust Tracker. (B) Particles are injected at different quantities into the test box and the count is noted 5 minutes after injection.

4.3 Experimental Rooms

Three different rooms were used in the study:

Simulation Room 1: This experimental was carried in a room of 17' x 13.6' x 9' designed to replicate an actual operation room equipped with a total of four or six sensors distributed as shown in **Figure 5:** i) two sensors MO-1 and DL-1 at 3 feet distance, ii) two sensors MO-2 and DL-2 or

one MO-2 at 6 feet and iii) two sensors MO-3 and DL-3 or DL-3 at 13 feet from the tested subject. The room had an air exchange rate of $\sim 20 \text{ h}^{-1}$. The measurements were taken every 1 minute. For MO sensors an integration time of 1 minute was used, for DL sensors an instant reading was registered every 1 minute. For each experimental a period of baseline was registered before starting the nebulization of 3 ml saline physiological solution, which was delivered through Piezoelectric Based Nebulizer: Pro-X Controller by Aerogen. Study cases I, II, III and V, which are typical healthcare facilities' applications were tested in this room. The specifications of instruments used in these cases is indicated below.

Simulation Room 2: This experimental was carried in a room 17' x 13.6' x 9' with the same particle sensor distribution of simulation room 1. The room had an air exchange rate of $\sim 2-3 \text{ h}^{-1}$, which is typically representative of a home environment. Case IV including the nebulization of 3 ml saline physiological solution with a mechanical pump nebulizer: VH Complete Compressor by Veridian under normal breathing conditions was tested in this room. The measurements were taken every 1 minute. For MO sensors and integration time of 1 minute was used, for DL sensors an instant reading was registered every 1 minute. For each experimental a period of baseline was registered before starting the nebulization.




Figure 5. Experimental Testing Plan to Assess the Efficacy.

Patient Room: This experimental was carried in a room of 17' x 13.6' x 9' equipped with a small bathroom of 4' x 6.25' x 9' in the back and an air ventilation rate, which is purposely set at 20 hour⁻¹ to minimize transmission of the disease via aerosols. A total of three sensors was distributed as shown in **Figure 1**: one DL-1 sensor at 3 feet distance, one DL-2 sensor at 6 feet and one DL-3 sensor at 9 feet from the COVID-19 patient. The room had an air exchange rate of ~20 h⁻¹. The 1-minute average reading were automatically recorded in the SD card of the instrument in the patient environment without disturbing the patient or putting the researchers at risk.

4.4 Real and Simulated Human Subjects

A total of 6 human test subjects participated in this study. All methods discussed comply with the proper regulations for proper reporting experimental on human subjects through the Office for Human Research Protections with approved assurance of compliance via Federal Wide Assurance (FWA number: ASU FWA 00009192). The human subjects were consented via ASU IRB: STUDY00006544 with the following IRB reviewers: Debra Murphy, empowered official and institutional official for IRB and Susan Metozky, IRB Compliance Office. In addition, we obtained informed consent from the subjects reported in the study for both study participation and publication in an online access publication. The human subjects were 22 to 55 years old, 5'5" to 6' tall, and had body mass indexes < 25. They participated in all of the tests in the study, except in tests corresponding to Case V, where a mannequin was used. The mannequin had an artificial lung simulator and ventilation was set to achieve 6 to 8 L per minute with a breath frequency of 16 breaths per minute, Tidal Volume (V_t) of 400 to 1000 mL.

4.5 Oxygen Delivery Methods

The oxygen delivery methods practiced in this study included the following specifications:

BiPAP: Full-Face Mask was used to deliver 21% of O_2 through non-invasive bilevel positive ventilation (NIPPV) ventilator set at inspiratory airway pressure (IPAP) of 12 cm of H_2O and Expiratory airway pressure (EPAP) of 5 cm of H_2O to test subjects. Resironics V60 by Philips provided a non-invasive bilevel positive airway pressure breathing support. Amara Full-Face Mask by Philips was used.

HFNC: Optiflow+ High flow nasal cannula (HFNC) from Fisher & Paykel was employed delivering 21% of O_2 and 30 or 60 L/min was used.

4.6 Aerosol Mitigation Methods

To mitigate the aerosol spread from the test subject the four different types of mitigations had the following specifications:

Biofilter: ThermoFlo filter by Arc Medical or Hudson Bacterial/viral Filter were employed, as it provides 99.9% bacterial and viral efficiency.

Disposable Surgical Mask: Common surgical mask used by doctors and patients.

Modified Silicone Disposable Mask: A custom designed mask by altering the conventional mask by Breezing™ (Tempe, AZ). The mask is fitted with the high flow nasal cannula or the mechanical nebulizer outlet port by slicing the silicon paddings of Breezing™ mask. The mask is also integrated with a biofilter in the exhalation tube followed by a battery-operated fan (Mouser Electronics).


Box Mitigation: A plastic box drilled two holes with each at opposite ends fitted with one-way air flow valve to provide one directional ventilation and this box is fitted with 2 MIL plastic curtains that runs throughout the test subject's body to reduce the leak of aerosols.

CHAPTER 5

MITIGATION SYSTEM EVALUATION

5.1 Validation test of the aerosol sensing instruments

Before starting the testing with human subjects, and the different oxygen delivery and mitigation methods, three MO sensors and one DL sensor were positioned in the simulation room at 3, 6, and 13 feet from a piezoelectric nebulizer with 3 mL of 0.19 μm polystyrene particles solution in ethanol, and a nebulization was delivered to the mannequin with a non-invasive positive pressure ventilator (NIPPV) and no aerosol mitigation system connected to it. **Figure 6** show the corresponding particle concentration profiles, indicating a clear peak after the initiation of the nebulizer, which corroborated the capacity of MO and DL sensors to capture aerosol plumes from nebulization within a room (simulation room 1) with a high air exchange rate (20 h^{-1}).






Figure 6. 0.19 μ m Physiological Solution Particle Count - concentration after nebulization of 3 ml saline at i. 3 feet, ii. 6 feet and iii. 13 feet from the nebulizer. (A) Particle count concentration levels from DY-01 placed at 3 feet from the nebulizer. (B) Particle count concentration levels at i) 3 feet (—), ii) 6 feet (---) and iii) 13 feet (----) from the nebulizer

5.2 Assessment of air exchange rate of the rooms


The air flowing in and out of the room is an important parameter as it is a major factor affecting the dispersion of aerosol. Increasing air exchange rates improves the room's capability of eliminating aerosols and can be determined by the declining slope of particles being vented out from the room¹⁵. To determine the air exchange rate in the tested room, the room was filled with


CO₂ gas, and the CO₂ gas concentration was monitored over time with a carbon dioxide sensor (Telaire by GE). The declining slope of recorded CO₂ levels yields air exchange rates of 20 hr⁻¹ for the simulation room 1 and patient room, and 2-3 hr⁻¹ for simulation room 2.

5.3 Description of study cases


The application of the mitigation methods to different oxygen therapies rendered five study cases described below.

Case I: The test subject used of a Full-Face Mask in bilevel positive airway pressure device (BPAP) therapy. The system was tested with and without the use of a biofilter located in its outlet (**Figure 3-I**). The nebulization therapy of 3 ml saline physiological solution was delivered through the Piezoelectric Based Nebulizer. The experimental were carried out in the simulation room 1 equipped with a total of four sensors: i) two sensors MO-1 and DL-1 at 3 feet distance, ii) one sensor MO-2 placed at 6 feet, and iii) one sensor MO-3 at 13 feet from the test subject. **Figure 7-I to 7-III** show detailed particle concentration profiles over time assessed at baseline with no nebulization, with nebulization and no mitigation, and with nebulization and mitigation at different distances. **Figure 7-IV** shows an overlapped particle concentration profiles over time assessed during nebulization and under no mitigation conditions, and a higher exposure to nebulizer aerosols are observed at 6 feet with respect to 3 feet location due to the air circulation streams in the room.





III



IV Full Face Mask with no Filter







Figure 7. Particle Count of BiPAP Therapy, with Full-Face Mask at (I) 3 feet, (II) 6 feet, and (III) 13 feet from the subject. The subject wore a Full Face Mask (FFM) with: i) no filter (—, —) , ii) filter and a leak of 30-60 L/min (---, ---) and iii) with filter and no leak (----, ----) from the mask. (IV) Comparative particle count profiles for Full Face Mask with no filter at different positions.


Case II: The test subject used of a High Flow Nasal Cannula (HFNC) therapy, with and without a surgical mask as mitigation system (**Figure 3-II**). The nebulization therapy of 3 ml saline physiological solution was delivered through the Piezoelectric Based Nebulizer. The experimental were carried out in the simulation room 1 equipped with a total of four sensors: i) two sensors MO-1 and DL-1 at 3 feet distance, ii) one sensor MO-2 placed at 6 feet, and iii) one sensor MO-3 at 13 feet from the test subject. **Figure 8-I to 8-III.** show detailed particle concentration profiles over time

assessed at baseline with no nebulization, with nebulization and no mitigation, and with nebulization and mitigation at different distances. **Figure 8-IV to 8-V** shows overlapped particle concentration profiles over time for all sensor locations without and with modified silicone mask + biofilter + fan use conditions, respectively. Similarly, to case I, a higher exposure to nebulizer aerosols is observed at 6 feet with respect to 3 feet location for the mask mitigation condition due to the air circulation streams in the room.















Figure 8. Particle Count of HFNC Therapy, with and without Surgical Mask at (I) 3 feet, (II) 6 feet, and (III) 13 feet from the subject. The subject wore: i) no mask (— , —) and ii) surgical mask (---, ---). (D) Comparative particle count profiles for no mask at different positions (IV) Comparative particle count profiles for no mask at different positions. (V) Comparative particle count profiles for the modified silicone mask at different positions.


Case III: The test subject used of a High Flow Nasal Cannula (HFNC) therapy, with and without a modified silicone mask and a biofilter and fan as mitigation system (**Fig. 3-III**). The nebulization therapy of 3 ml saline physiological solution was delivered through the Piezoelectric Based Nebulizer. The experimental were carried out in the simulation room 1 equipped with a total of six sensors: i) two sensors MO-1 and DL-1 at 3 feet distance, ii) one sensor MO-2 placed at 6 feet, and iii) one sensor MO-3 at 13 feet from the test subject. **Figure 9-I to 9-III** show detailed particle

concentration profiles over time assessed at baseline with no nebulization, with nebulization and no mitigation, and with nebulization and mitigation at different distances. **Figure 9-IV to 9-V** shows an overlapped particle concentration profiles over time for all sensor locations without and with modified silicone mask + biofilter + fan use conditions, respectively. **Figure 9-V** (case of use of the modified mask+ biofilter + fan) shows the clear reduction of the nebulizer aerosol particles close to the baseline levels.















Figure 9. Particle Count of HFNC Therapy, with and without Modified Silicone Mask at (A) 3 feet, (B) 6 feet, and (C) 13 feet from the subject. The subject wore: i) no mask (—, —) and ii) modified silicone mask (---, ---). (D) Comparative particle count profiles for no mask at different positions. (E) Comparative particle count profiles for the modified silicone mask at different positions.


Case IV: The test subject was breathing naturally without the help of any oxygen therapy and was tested with and without the modified silicone mask + biofilter + fan as mitigation system (**Fig. 3-IV**). The nebulization therapy of 3 ml saline physiological solution was delivered through the Mechanical Based Nebulizer. This nebulizer delivered a 4L/min air with the aerosol. The experimental were carried out in the simulation room 2 equipped with a total of six sensors: i) two sensors MO-1 and

DL-1 at 3 feet distance, ii) two sensors MO-2 and DL-2 placed at 6 feet, and iii) two sensors MO-3 and DL-3 at 13 feet from the test subject. **Figure 10-I to 10-III** show detailed particle concentration profiles over time assessed at baseline with no nebulization, with nebulization and no mitigation, and with nebulization and mitigation at different distances. **Figure 10-IV to 10-V** shows an overlapped particle concentration profiles over time for all sensor locations without and with modified silicone mask + biofilter + fan use conditions, respectively. **Figure 10-V** (case of use of the modified mask+ biofilter + fan.) shows the clear reduction of the nebulizer aerosol particles close to the baseline levels.














Figure 10. Particle Count of Natural Breathing, with and without Modified Silicon Mask at (I) 3 feet, (II) 6 feet, and (III) 13 feet from the subject. The subject wore: i) no mask (— , —) and ii) modified silicone mask (---, ---). (IV) Comparative particle count profiles for no mask at different positions. (V) Comparative particle count profiles for the modified silicone mask at different positions.

Case V: The mannequin was ventilated with a Full-Face Mask connected to a Non-Invasive Positive Pressure Ventilator (NIPPV) and was tested with and without a mitigation box that had a biofilter attached to it as mitigation system (**Fig. 3-V**). The mitigation box design was inspired on the isolation hood that has gained attention with the focus of mitigation SARS-CoV-2 pathogens²³. The nebulization therapy of 3 ml saline physiological solution was delivered through the


Piezoelectric Based Nebulizer. The experimental were carried out in the simulation room 1 equipped with a total of six sensors: i) two sensors MO-1 and DL-1 at 3 feet distance, ii) two sensors MO-2 and DL-2 placed at 6 feet, and iii) two sensors MO-3 and DL-3 at 13 feet from the test subject. **Figure 11-I to 11-III** show detailed particle concentration profiles over time assessed at baseline with no nebulization, with nebulization and no mitigation, and with nebulization and mitigation at different distances. **Figure 11-IV to 11-V** shows an overlapped particle concentration profiles over time for all sensor locations with and without modified mitigation box, respectively.






III

13 Feet








Figure 11. Particle Count of BiPAP Therapy, with Mitigation Box at (I) 3 feet, (II) 6 feet, and (III) 13 feet from the subject. The subject wore: i) Box mitigation (—, —) and ii) No mask (---, - - -). (IV) Comparative particle count profiles for box mitigation at different positions. (V) Comparative particle count profiles for no mask at different positions.


CHAPTER 6

RESULTS

Figure 12-I to 12-II shows a summary of the efficacy of different mitigation systems that are applied to different oxygen therapies and nebulizers. The efficacy of the nebulization therapy mitigation is assessed using the peak and the area under the aerosol particle counts / feet³ transient curves for both the no mitigation and mitigation cases. As an example, **Figure 12-III** shows the aerosol particle counts / feet³ transient over time at 3 feet from a patient during a nebulization therapy with and without the mitigation system. The following equations are used in connection with transients as the ones shown in **Figure 12-III** for calculation of the peak and area efficacy of the systems:

$$\text{Peak Efficacy (\%)} = \frac{[(\text{Peak with no mitigation solution}) - (\text{Peak with mitigation solution})]}{\text{Peak with no mitigation solution}} \quad (\text{Eq. 1})$$

$$\text{Area Efficacy (\%)} = \frac{[(\text{Area with no mitigation solution}) - (\text{Area with mitigation solution})]}{\text{Area with no mitigation solution}} \quad (\text{Eq. 2})$$



II




Figure 12. Efficacy of Aerosol Mitigation methods applied in this work calculated from (I) Peak of particle count concentration (II) Area under the curve of particle count concentration vs. time during the simulated nebulization therapies. Reference: (**Case I**): BPAP therapy with Full-Face Mask, (**Case II**) HFNC therapy with surgical mask, (**Case III**) HFNC therapy with modified silicone mask, (**Case IV**) Natural breathing with modified silicone mask, (**Case V**) NIPPV therapy with box mitigation

The peak and area efficacies are evaluated for curves taken at 3, 6 and 13 feet with experimental setup similar to the one shown in **Figure 1** and averaged to report the efficacy registered for the particular condition. **Figure 7-11** show the corresponding particle counts / feet³ profiles obtain for *Case I* (**Figure 7**), *Case II* (**Figure 8**), *Case III* (**Figure 9**), *Case IV* (**Figure 10**),

and Case V (**Figure 11**). As it can be observed in **Figure 12-I to 12-II**, Case II and Case V, which includes a surgical mask and a mitigation box, have only a mitigation efficacy of 10-24% to capture of nebulizer aerosols. Unfortunately, these have been commonly used methods to mitigate nebulization's aerosol dispersions. On the contrary, Cases I, III and IV render an efficacy of 98-100% to capture of nebulizer aerosols. Since HFNC and BPAP are the preferred oxygen therapies to delay need for invasive mechanical ventilation in COVID-19 patients, Cases I and III includes representative solutions to mitigated aerosols during nebulization of COVID-19 patients or any patient with an infectious respiratory disease. Further, in order to produce 100% efficacy in Case I, it is observed that the complete sealing of the Full-Face Mask is required, given that any leak generates aerosol release (**Figure 7-I to 7-III**). Cases III and IV include the use of a modified silicone disposable mask which has: 1- an inlet port for the delivery of the oxygen delivery therapy, 2- an outlet port for eliminating the residual aerosols from the therapy together with the unused air from the therapy, which is equipped with the filtering membrane and the fan, and 3- a safety/emergency inlet port with a one-way valve for inhalation of ambient air to be used in case the oxygen therapy stops. The mask incorporates the synergic action of convective transport from a fan and the filtration process from a biofilter with bacterial and viral filtering membrane of 0.2 μm pore size and multiple layers of microfiber material (**Figure 13**). In addition, it features an ergonomically designed silicone edge that enables adaptation to the user's face, which is essential for a good seal to prevents leakage of contaminated air. The fan allows the driving of the aerosols through the biofilter, assuring all aerosols of 0.2 μm and larger are captured before the main air stream of the therapy flow exists the system to the atmospheric environment. It is driven by a 5V portable and rechargeable battery that can sustain the fan's driving under the system's planned operative conditions for ~6 hours continuously. Considering a nebulization therapy (including nebulization + therapy) takes 20 minutes, the battery can sustain 12 therapies.

Mask* - Nebulizer Aerosol Mitigation System




Figure 13. Modified Silicone Mask

CHAPTER 7

CONCLUSION

In summary, we demonstrate here existing masks and boxes mitigation methods are inefficient for capturing aerosols generated from nebulizers. On the contrary, we probed the high efficacy for nebulization aerosol capture of a biofilter used in conjunction with a fully sealed Full-Face Mask (FFM) in **BiPAP** therapy and the modified silicone disposable mask in **HFNC** therapy and **Normal Breathing**. Given the crisis the whole world is facing with SARS-CoV-2 spread, the highly efficient aerosol trapping systems can help to limit transmission of the airborne virus in the current COVID-19 epidemic or other current respiratory infectious disease pathogens or future pandemics.

REFERENCES

1. Galbadage, T., Peterson, B.M. & Gunasekera, R.S. Does COVID-19 spread through droplets alone? *Front. Public Health* 24 April, <https://doi.org/10.3389/fpubh.2020.00163> (2020).
2. Zayas, G. et al. Cough aerosol in healthy participants: fundamental knowledge to optimize droplet-spread infectious respiratory disease management. *Bmc Pulm Med* 12 (2012).
3. Santarpia, J.L. et al. Aerosol and Surface Transmission Potential of SARS-CoV-2. *medRxiv* 2020.03.23., <https://doi.org/10.1101/2020.1103.1123.20039446> (2020).
4. Sunjaya, A.P. & Jenkins, C. Rationale for universal face masks in public against COVID-19. *Respirology* 25, 678-679 (2020).
5. Li, Y. et al. Role of ventilation in airborne transmission of infectious agents in the built environment – a multidisciplinary systematic review. *Indoor Air, International Journal of Indoor Environment and Health* 17, 2-18 (2007).
6. Escombe, A.R., Ticona, E., Chávez-Pérez, V., Espinoza, M. & Moore, D.A.J. Improving natural ventilation in hospital waiting and consulting rooms to reduce nosocomial tuberculosis transmission risk in a low resource setting. *BMC Infect Dis.* 19, 88. Published 2019 Jan 2025. doi:2010.1186/s12879-12019-13717-12879 (2019).
7. Correia, G., Rodrigues, L., Gameiro da Silva, M. & Gonçalves, T. Airborne route and bad use of ventilation systems as non-negligible factors in SARS-CoV-2 transmission. *Med Hypotheses.* 141, doi:10.1016/j.mehy.2020.109781 (2020).
8. Tran, K., Cimon, K., Severn, M., Pessoa-Silva, C.L. & Conly, J. Aerosol generating procedures and risk of transmission of acute respiratory infections to healthcare workers: a systematic review. *PLoS One*, <https://doi.org/10.1371/journal.pone.0035797> (2012).
9. Schunemann, H.J. et al. Ventilation Techniques and Risk for Transmission of Coronavirus Disease, Including COVID-19. *Annals of Internal Medicine* 173 (2020).
10. WHO WHO International Weekly Epidemiological Update. World Health Organization 3 November 2020. Accessed November 7, 2020. <https://www.who.int/publications/m/item/weekly-epidemiological-update---3-november-2020> (2020).
11. I, A. & Ysrafil Severe Acute Respiratory Syndrome Coronavirus 2 (SARS-CoV-2): An overview of viral structure and host response. *Diabetes Metab Syndr.* 14, 407-412 (2020).
12. CDC CDC COVID Data Tracker: Maps, charts, and data provided by CDC. Coronavirus Disease 2019 (COVID-19). https://covid.cdc.gov/covid-data-tracker/?CDC_AA_refVal=https%3A%2F%2Fwww.cdc.gov%2Fcoronavirus%2F2019-ncov%2Fcases-updates%2Fcases-in-us.html#cases_casesper100klast7days (2020)

13. Vuorinen, V. et al. Modelling aerosol transport and virus exposure with numerical simulations in relation to SARS-CoV-2 transmission by inhalation indoors. *Safety Science* 130, 104866 (2020).
14. Liu, Y. et al. Aerodynamic analysis of SARS-CoV-2 in two Wuhan hospitals. *Nature* (2020).
15. Kang, S. et al. Crystal structure of SARS-CoV-2 nucleocapsid protein RNA binding domain reveals potential unique drug targeting sites. *Acta Pharmaceutica Sinica B* 10, 1228-1238 (2020).
16. Fennelly, K.P. Particle sizes of infectious aerosols: implications for infection control. *Lancet Respiratory Medicine* 8 (2020).
17. Baron, P. Generation and behavior of airborne particles. https://www.cdc.gov/niosh/topics/aerosols/pdfs/aerosol_101.pdf
18. Stadnytskyi, V., Bax, C.E., Bax, A. & Anfinrud, P. The airborne lifetime of small speech droplets and their potential importance in SARS-CoV-2 transmission. *P Natl Acad Sci USA* 117, 11875-11877 (2020).
19. Mousavi, E.S., Kananizadeh, N., Martinello, R.A. & Sherman, J.D. COVID-19 Outbreak and Hospital Air Quality: A Systematic Review of Evidence on Air Filtration and Recirculation. *Environ. Sci. Technol.* (2020).
20. CDC Scientific Brief: SARS-CoV-2 and Potential Airborne Transmission. *Coronavirus Disease 2019 (COVID-19)* Oct. 5th (2020).
21. Simones, M.P. et al. Measurement of the size and charge distribution of sodium chloride particles generated by an Aerobe Pro® pharmaceutical nebulizer. *European Journal of Nanomedicine* Volume 6: Issue 1 (22 Jan 2014).
22. Moix, A.L.S. et al. Physicochemical Aspects and Efficiency of Albuterol Nebulization: Comparison of Three Aerosol Types in an In Vitro Pediatric Model. *RESPIRATORY CARE* 60, 38-46 (2015).
23. Shaw, K.M. et al. Intensive care unit isolation hood decreases risk of aerosolization during noninvasive ventilation with COVID-19. *Can J Anesth* (2020).